

MAXI+3 WHEEL STAND WITH SET OF 3 CUFFS (S), (M), (L)

Ref. 522708



CE
0123



MAXI+3 WHEEL STAND WITH SET OF 3 CUFFS (S), (M), (L)



Ref. 522708

SPENGLER technology in a blood pressure monitor with a dial with a diameter of 147mm.

Description

HIGH MEASUREMENT ACCURACY

- Dual hoses for linear pressure build-up

COMFORT

- Ambidextrous
- Remote-readable display
- 2.5-meter coiled hose
- Adjustable from 95 cm to 135 cm
- Ergonomic, large-volume bulb for rapid inflation
- Storage basket behind the dial (available as an option)
- Comes with a wheeled stand adjustable from 95 cm to 135 cm

CUFF

- Comes with a set of 3 cuffs (S, M, L), nylon

MAXI+3 WHEEL STAND WITH SET OF 3 CUFFS (S), (M), (L)

Ref. 522708



Features

CE Marking	CE0123	Medical device	Class Im (measuring function)
Bar code	3700446031169	HS Code	90189010
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	4.36
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		