

**S/S BOWL WITHOUT SPOUT, FLAT
BOTTOM, DIAM 100 MM, 300 ML, WITH
COVER**

Ref. AI99101



19/06/2026 10:03:42 – For information purposes only, can be subject to change

S/S BOWL WITHOUT SPOUT, FLAT BOTTOM, DIAM 100 MM, 300 ML, WITH COVER



Ref. [AI99101](#)

Allows surgical instruments to be stored and transported

Description

Holtex instruments are all made of premium stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to ensure good durability and corrosion resistance.

The steels come into contact with human tissues, so they must comply with international standards and regulations.

Holtex guarantees you medical devices that comply with Directive 93/42/EC and ISO 13485 (2016).

- In the "ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" leaflet

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104740495	HS Code	90189084
Warranty (months)	12	Country	Pakistan

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	21	Bare Unit / Width (cm)	21.5
Packed unit / Height (cm)	4	Bare Unit / Height (cm)	4.6
Packed unit / Length (cm)	38	Bare Unit / Length (cm)	38
Packed unit / Weight	0.1	Bare Unit / Weight (kg)	0.096

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---