

CHAPUT FORCEPS 16 CM

Ref. IP11616




CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
2 ANS
2 YEARS
GUARANTEE

STAINLESS
INOX
STEEL



19/06/2026 03:50:52 – For information purposes only, can be subject to change

CHAPUT FORCEPS 16 CM

Ref. IP11616



Features

CE Marking	CE	Medical device	Class I
Bar code	08961104716186	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.028
Bare Unit / Weight (kg)	0.028		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 19/06/2026 03:50:52 – For information purposes only, can be subject to change