

DOYEN ABDOMINAL RETRACTOR - 24 CM - 35 MM X 85 MM



Ref. IV03486

Features

CE Marking	CE	Medical device	Class I
Bar code	08961104722323	HS Code	90189084
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.234
Bare Unit / Weight (kg)	0.234		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 19/06/2026 07:57:39 – For information purposes only, can be subject to change