

ELECTRODES FOR TRAINING DEFIBRILLATOR (5 PAIRS)

Ref. SJACCE69



Description

Designed solely for training in the use of a defibrillator in cardiac resuscitation. Do not use on a human.

Features

CE Marking	NONE	Medical device	Class I
Bar code	3661040137285	HS Code	90181990
Warranty (months)	12	Country	Italy

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	16
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	22
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	23
Packed unit / Weight	0.1	Bare Unit / Weight (kg)	0.1

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---