

MAGISTER STETHOSCOPE FOR CHILDREN / BLACK

Ref. 507210



22/06/2026 22:56:15 – For information purposes only, can be subject to change

MAGISTER STETHOSCOPE FOR CHILDREN / BLACK

Ref. 507210



=====

Description

For pediatric auscultation. Stainless steel dual head chestpiece for pediatric use. Multifrequency diaphragm. Angled binaural with built-in spring and flexible, screw-on tips. Equipped with a non-chill ring for little ones' comfort.

Available in :

Blue : Ref 507214

Black: Ref 507210

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009120	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	4
Packed unit / Length (cm)	32	Packed unit / Weight	0.375
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-10	Max temperature	40
Minimum humidity	NA	Maximum humidity	95
OTHER FEATURES			
Colors	Black	Packed unit / Quantity	1