

PEAK FLOWMETER

Ref. L36147



CE
0197



Peak flow meter designed for accurate monitoring of respiratory function.

Description

Un débitmètre offert pour 5 boîtes d'embouts

22/06/2026 23:48:28 – For information purposes only, can be subject to change

PEAK FLOWMETER

Ref. L36147



Features

CE Marking	CE0197	Medical device	Class Im (measuring function)
Bar code	3700446001780	HS Code	90261081
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	38	Packed unit / Height (cm)	38
Packed unit / Length (cm)	50	Packed unit / Weight	0.075
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	1	Max temperature	50
Minimum humidity	10	Maximum humidity	95
OTHER FEATURES			
Code LPP	6115365	Packed unit / Quantity	1