

=====

Ref. BCTOCO33

COMeN



CE
1639



=====

Description

=====

 24/06/2026 01:40:03 – For information purposes only, can be subject to change

=====

COMeN

Features

CE Marking	CE1639	Medical device	Class IIa
Bar code	6945429082487	HS Code	9018908400
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES