

=====

Ref. BCTOCO36

COMEN



CE  
1639



=====

## Description

=====

📅 24/06/2026 00:59:18 – For information purposes only, can be subject to change

=====

**COMeN**

## Features

|                   |               |                |            |
|-------------------|---------------|----------------|------------|
| CE Marking        | CE1639        | Medical device | Class I    |
| Bar code          | 6945429083880 | HS Code        | 9018908400 |
| Warranty (months) | 6             | Country        | China      |

### DIMENSIONS AND WEIGHT

|                           |    |                           |    |
|---------------------------|----|---------------------------|----|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | nc |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

### OTHER FEATURES