

Ref. 332012



CE



26/06/2026 00:12:57 – For information purposes only, can be subject to change

Ref. 332012


Description

Cotton cuff with highly resistant rubber bladder.
 Twin 7,5 cm tubing fitted with quick CPC coupling. Legal markings printed.
 Washable at 30-40°C.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446028480	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.21
Bare Unit / Weight (kg)	0.21		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 26/06/2026 00:12:57 – For information purposes only, can be subject to change