

Ref. 332029



CE



25/06/2026 23:50:17 – For information purposes only, can be subject to change

Ref. 332029


Description

Cotton cuff with highly resistant rubber bladder.
 Twin 7,5 cm tubing fitted with quick CPC coupling. Legal markings printed.
 Washable at 30-40°C.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446008796	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.094
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 25/06/2026 23:50:17 – For information purposes only, can be subject to change