

**Ref. 332039**



CE



26/06/2026 00:06:59 – For information purposes only, can be subject to change

---

**Ref. 332039**


## Description

Cotton cuff with highly resistant rubber bladder.  
 Twin 7,5 cm tubing fitted with quick CPC coupling. Legal markings printed.  
 Washable at 30-40°C.

## Features

|                   |               |                |            |
|-------------------|---------------|----------------|------------|
| CE Marking        | CE            | Medical device | Class I    |
| Bar code          | 3700446008789 | HS Code        | 9018901000 |
| Warranty (months) | 6             | Country        | Indonesia  |

### DIMENSIONS AND WEIGHT

|                           |    |                           |       |
|---------------------------|----|---------------------------|-------|
| Packed unit / Width (cm)  | 15 | Packed unit / Height (cm) | 2     |
| Packed unit / Length (cm) | 13 | Packed unit / Weight      | 0.062 |

### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|

 26/06/2026 00:06:59 – For information purposes only, can be subject to change