

MAGISTER STETHOSCOPE FOR ADULT / BLACK

Ref. 507202



CE



28/06/2026 02:14:47 – For information purposes only, can be subject to change

MAGISTER STETHOSCOPE FOR ADULT / BLACK

Ref. 507202



Description

For general-purpose auscultation. Stainless steel dual head chestpiece for adults. Multifrequency diaphragm. Angled binaural with built-in spring and flexible, screw-on tips. Ideal grip and high performance make it the perfect product for general and specialist use.

Available in en Grey (ref 507204), Red (ref: 507201), Green (ref :507209), Bleu (ref : 507202), Black (ref: 507200)

STORAGE CONDITIONS :

Atmospheric pressure : 700-1060 hPa

Temperature: -10° to 40° C

Humidity: 95% max

MAGISTER STETHOSCOPE FOR ADULT / BLACK



Ref. 507202

Features

| | | | |
|------------------------------------|---------------|---------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3700446009076 | HS Code | 90181990 |
| Warranty (months) | 60 | | |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 16 | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | 32 | Packed unit / Weight | 0.398 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | -10 | Max temperature | 40 |
| Minimum humidity | 95 | Maximum humidity | 95 |
| OTHER FEATURES | | | |
| Colors | Navy blue | Packed unit / Quantity | 1 |

28/06/2026 02:14:47 – For information purposes only, can be subject to change