


# MAGISTER STETHOSCOPE FOR ADULT / BLACK

Ref. 507209



CE



28/06/2026 01:41:07 – For information purposes only, can be subject to change

# MAGISTER STETHOSCOPE FOR ADULT / BLACK

Ref. 507209



## Description

For general-purpose auscultation. Stainless steel dual head chestpiece for adults. Multifrequency diaphragm. Angled binaural with built-in spring and flexible, screw-on tips. Ideal grip and high performance make it the perfect product for general and specialist use.

Available in en Grey ( ref 507204), Red (ref: 507201), Green (ref :507209), Bleu (ref : 507202), Black (ref: 507200)

### STORAGE CONDITIONS :

Atmospheric pressure : 700-1060 hPa

Temperature: -10° to 40° C

Humidity: 95% max

# MAGISTER STETHOSCOPE FOR ADULT / BLACK



Ref. 507209

## Features

CE Marking	CE	Medical device	Class I
Bar code	3700446004002	HS Code	9018.90.20.00-9
Warranty (months)	60		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.398
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-10	Max temperature	40
Minimum humidity	95	Maximum humidity	95
<b>OTHER FEATURES</b>			
Colors	Green	Packed unit / Quantity	1

28/06/2026 01:41:07 – For information purposes only, can be subject to change