


# POLYESTER AUTOCLAVABLE CUFF STRAIGHT WITH 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GBS101

**DDM** 



CE



📅 01/07/2026 01:27:41 – For information purposes only, can be subject to change

# POLYESTER AUTOCLAVABLE CUFF STRAIGHT WITH 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET

Ref. GBS101



Those cuffs are autoclavable!

## Description

Cuffs for DDM pneumatic tourniquets are made of polyester. They are disinfectable and autoclavable and composed of a removable silicon bladder for easy autoclaving.

Those ones are designed in straight form

Several sizes are available

### Dimensions of the bladder in cm

|                           |             |
|---------------------------|-------------|
| Cuff infant (GBS101)      | 19 x 4 cm   |
| Cuff child (GBS102)       | 26.5 x 5 cm |
| Cuff small adult (GBS101) | 35.5 x 6 cm |
| Cuff adult (GBS103)       | 53 x 6 cm   |

# POLYESTER AUTOCLAVABLE CUFF STRAIGHT WITH 1 BLADDER FOR DDM PNEUMATIC TOURNIQUET



Ref. GBS101

## Features

|                   |               |                |          |
|-------------------|---------------|----------------|----------|
| CE Marking        | CE            | Medical device | Class I  |
| Bar code          | 3700468203476 | HS Code        | 90189084 |
| Warranty (months) | 24            | Country        | France   |

### DIMENSIONS AND WEIGHT

|                           |    |                           |    |
|---------------------------|----|---------------------------|----|
| Packed unit / Width (cm)  | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight      | nc |

### SPECIFIC STORAGE CONDITIONS

|                     |     |                  |    |
|---------------------|-----|------------------|----|
| Minimum temperature | -10 | Max temperature  | 40 |
| Minimum humidity    | 30  | Maximum humidity | 40 |

### OTHER FEATURES