

## MAGISTER STETHOSCOPE FOR ADULT / BLACK

Ref. 507200



08/07/2026 00:47:16 – For information purposes only, can be subject to change

## MAGISTER STETHOSCOPE FOR ADULT / BLACK

Ref. 507200



=====

### Description

For general-purpose auscultation. Stainless steel dual head chestpiece for adults. Multifrequency diaphragm. Angled binaural with built-in spring and flexible, screw-on tips. Ideal grip and high performance make it the perfect product for general and specialist use.

Available in en Grey ( ref 507204), Red (ref: 507201), Green (ref :507209), Bleu (ref : 507202), Black (ref: 507200)

#### STORAGE CONDITIONS :

Atmospheric pressure : 700-1060 hPa

Temperature: -10° to 40° C

Humidity: 95% max

## MAGISTER STETHOSCOPE FOR ADULT / BLACK

Ref. 507200



### Features

|                             |               |                           |          |
|-----------------------------|---------------|---------------------------|----------|
| CE Marking                  | CE            | Medical device            | No       |
| Bar code                    | 3700446009069 | HS Code                   | 90181990 |
| Warranty (months)           | 60            |                           |          |
| DIMENSIONS AND WEIGHT       |               |                           |          |
| Packed unit / Width (cm)    | 16            | Packed unit / Height (cm) | 4        |
| Packed unit / Length (cm)   | 32            | Packed unit / Weight      | 0.398    |
| SPECIFIC STORAGE CONDITIONS |               |                           |          |
| Minimum temperature         | -10           | Max temperature           | 40       |
| Minimum humidity            | 95            | Maximum humidity          | 95       |
| OTHER FEATURES              |               |                           |          |
| Colors                      | Black         | Packed unit / Quantity    | 1        |