

**Ref. 336401**



CE



## Description

Nylon cuffs with quick coupling for Lian® with 60 cm single tubing bladder fitted with a quick coupling for Lian® NM, Lian® Métal, Lian® Nano, Lian® Scénic, Ergo Sphyg®, Lian® Color. Washable at 30-40° C.

09/07/2026 01:12:47 – For information purposes only, can be subject to change

-----  
**Ref. 336401**


## Features

CE Marking	CE	Medical device	Class I
Bar code	3700446026615	HS Code	9018901000
Warranty (months)	6	Country	Indonesia
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.29
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 09/07/2026 01:12:47 – For information purposes only, can be subject to change