

**MEDIPLUGS OPERATING ROOM SHOES,  
NON-PERFORATED, 35, NAVY BLUE**

Ref. WSAM\_43BC

medi **Plogs**



📅 09/07/2026 01:22:28 – For information purposes only, can be subject to change

## MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE

Ref. WSAM\_43BC



=====

### Description

Clogs suitable for modern hospitals, operating theatres and intensive care units.

Products designed and manufactured in collaboration with researchers and doctors.

- Antistatic clogs with non-slip and non-chafing soles
- The back of the heel is sloped and rounded to make walking easier
- Ventilation system allows air to circulate
- Toe prints to enhance hold and spurs to aid the circulation
- 100% CFC-free polyurethane
- Washable by hand or machine at 60°C
- Sizes: 35 - 47

## MEDIPLUGS OPERATING ROOM SHOES, NON-PERFORATED, 35, NAVY BLUE

Ref. WSAM\_43BC



### Features

|                                    |               |                           |          |
|------------------------------------|---------------|---------------------------|----------|
| CE Marking                         | NONE          | Medical device            | No       |
| Bar code                           | 3661040100173 | HS Code                   | 64029993 |
| Country                            | Germany       |                           |          |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |          |
| Packed unit / Width (cm)           | nc            | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)          | nc            | Packed unit / Weight      | 1.16     |
| Bare Unit / Weight (kg)            | 1             |                           |          |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |          |
| Minimum temperature                | NA            | Max temperature           | NA       |
| Minimum humidity                   | NA            | Maximum humidity          | NA       |
| <b>OTHER FEATURES</b>              |               |                           |          |
| Colors                             | Light blue    | Shoe Size                 | 43       |