


CARDIOMATE EVI 7" TABLET ECG SPENGLER

Ref. 260050



CE
0123



18/07/2026 00:02:42 – For information purposes only, can be subject to change

CARDIOMATE EVI 7" TABLET ECG SPENGLER

Ref. 260050



Spengler expands its range of ECGs and offers you its latest innovation: the CARDIOMATE EVI 7" tablet, a fully autonomous technological product.

Description

Very easy to use:

- 7" color touch screen
- Portable: compact size and light weight, easy to carry: 0.8kg and 197 x 112.4 x 26.1 mm
- Large battery capacity, rechargeable on mains
- Allows to print ECGs directly on a printer via wifi network

High-performance technology:

- High accuracy: sampling rate 32,000 Hz
- Detection and marking of the pacemaker
- Automatic interpretation (Glasgow) and calculation of intervals
- Input of patient information using a keyboard or barcode reader
- real time visualization of the 12 tracks at the same time as automatic interpretation and interval calculations. (possibility of a 18 leads tablet optional)
- Data export via wifi or USB port under .pdf or .xml

**CARDIOMATE EVI 7" TABLET ECG
SPENGLER****Ref. 260050****Features**

CE Marking	CE0123	Medical device	Class IIa
Bar code	6971049341789	HS Code	9018110000
Warranty (months)	36	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	10	Packed unit / Height (cm)	28
Packed unit / Length (cm)	33	Packed unit / Weight	2.8
Bare Unit / Weight (kg)	0.8		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	95

OTHER FEATURES