

REPLACEMENT INFRARED EAR THERMOMETER FOR NC3/NC5

Ref. BCTHER01

COMeN



The wireless infrared thermometer that gives a measurement in 2 to 3 seconds!

REPLACEMENT INFRARED EAR THERMOMETER FOR NC3/NC5

COMEN

Ref. BCTHER01

Features

CE Marking	NONE	Medical device	Class IIa
Bar code	3661040219349	HS Code	9025199010
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	17	Packed unit / Height (cm)	6
Packed unit / Length (cm)	22	Packed unit / Weight	nc
Bare Unit / Weight (kg)	0.1		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--