



## EMERAIR OXYGEN EMERGENCY BAG, RED

Ref. MEEME03\_R

**ELITE BAGS**



27/07/2026 23:50:50 – For information purposes only, can be subject to change

## EMERAIR OXYGEN EMERGENCY BAG, RED

Ref. MEEME03\_R



**For emergency services specialized in first aid**

### Description

- Interior totally adaptable to the needs of the user thanks to extractable velcro separators
- Elastic bands in all interior flat surfaces
- Reflective strips night and day- 2 external pockets of medium size, 2 external pockets of large capacity, with elastic bands
- 4 extractable compartments with the transparent front part differentiated according to the standardized color code (blue: respiratory system, red: circulatory system, green: various surgical equipment, yellow: pediatric equipment)
- Isothermal amp maker with 27 elastics
- Shoulder strap, handles that surround the bag for better resistance, concealed backpack straps and banana
- Material: 1000D polyamide
- Weight: 3.65 kg- Outdoor measures: 30 x 50 x 23 cm / 32.96 L
- Recommended max. capacity: 25 kg
- Max. capacity of an oxygen cylinder 3L
- Delivered with an ampoulier
- Optional trolley : MECHAR02
- Available in waterproof version
- Sold empty

# EMERAIR OXYGEN EMERGENCY BAG, RED



Ref. MEEME03\_R

## Features

|                                    |               |                         |            |
|------------------------------------|---------------|-------------------------|------------|
| CE Marking                         | NONE          | Medical device          | No         |
| Bar code                           | 8436564540133 | HS Code                 | 4202929890 |
| Warranty (months)                  | 24            | Country                 | China      |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                         |            |
| Packed unit / Width (cm)           | 32            | Bare Unit / Width (cm)  | 30         |
| Packed unit / Height (cm)          | 19            | Bare Unit / Height (cm) | 23         |
| Packed unit / Length (cm)          | 50            | Bare Unit / Length (cm) | 50         |
| Packed unit / Weight               | 4.1           | Bare Unit / Weight (kg) | 3.65       |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                         |            |
| Minimum temperature                | NA            | Max temperature         | NA         |
| Minimum humidity                   | NA            | Maximum humidity        | NA         |
| <b>OTHER FEATURES</b>              |               |                         |            |
| Packed unit / Quantity             | 1             |                         |            |