

Ref. C1.416.944

LUXAMED®



23/07/2026 22:34:05 – For information purposes only, can be subject to change

=====

Ref. C1.416.944

LUXAMED®

=====

Description

Available in several colours:

White (C1.416.914)

Light blue (C1.416.934)

Fuchsia (C1.416.974)

Yellow (C1.416.944)

Kids (C1.816.914)

Black (C1.416.114)

Green (C1.416.954)

Purple (C1.416.964)

Features

CE Marking	NONE	Medical device	No
Bar code	04260381963790	HS Code	90189084
Warranty (months)	24	Country	Germany

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90

OTHER FEATURES

Colors	YELLOW
--------	--------