

## DIAGNOSTIC PENLIGHT - WITH STANDARD BULB - LUXAMED

Ref. D1.212.212

**LUXAMED®**



23/07/2026 00:29:38 – For information purposes only, can be subject to change

## DIAGNOSTIC PENLIGHT - WITH STANDARD BULB - LUXAMED

**LUXAMED®**Ref. D1.212.212

Penlight for eye and throat diagnostic., long life technology.

### Description

- Top quality ABS plastic housing
- Warm white light : comfortable for eyes
- Operating voltage 3.0 V
- Incl. 2 AAA batteries

Available in :

- White: D1.212.912
- Blue: D1.212.212
- Grey: D1.212.312
- Yellow: D1.212.512
- Kids : D1.312.912
- Black: D D1.212.112
- Red : D1.212.412
- Green : D1.212.612

# **DIAGNOSTIC PENLIGHT - WITH STANDARD BULB - LUXAMED**

**LUXAMED®**
Ref. D1.212.212

## **Features**

|                                    |                |                           |          |
|------------------------------------|----------------|---------------------------|----------|
| CE Marking                         | CE             | Medical device            | Class I  |
| Bar code                           | 04260381961758 | HS Code                   | 90189084 |
| Warranty (months)                  | 24             | Country                   | Germany  |
| <b>DIMENSIONS AND WEIGHT</b>       |                |                           |          |
| Packed unit / Width (cm)           | nc             | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)          | nc             | Packed unit / Weight      | nc       |
| Bare Unit / Weight (kg)            | 0.04           |                           |          |
| <b>SPECIFIC STORAGE CONDITIONS</b> |                |                           |          |
| Minimum temperature                | -10            | Max temperature           | 40       |
| Minimum humidity                   | 10             | Maximum humidity          | 90       |
| <b>OTHER FEATURES</b>              |                |                           |          |
| Colors                             | Blue           |                           |          |