

=====

**Ref. E1.329.012**

**LUXAMED®**



CE



=====

## Description

=====

📅 28/07/2026 00:39:16 – For information purposes only, can be subject to change

=====

**Ref. E1.329.012**

**LUXAMED®**

## Features

|                                    |                |                           |          |
|------------------------------------|----------------|---------------------------|----------|
| CE Marking                         | CE             | Medical device            | Class I  |
| Bar code                           | 04260381960430 | HS Code                   | 90189084 |
| Warranty (months)                  | 24             | Country                   | Germany  |
| <b>DIMENSIONS AND WEIGHT</b>       |                |                           |          |
| Packed unit / Width (cm)           | nc             | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)          | nc             | Packed unit / Weight      | nc       |
| <b>SPECIFIC STORAGE CONDITIONS</b> |                |                           |          |
| Minimum temperature                | 5              | Max temperature           | 40       |
| Minimum humidity                   | 10             | Maximum humidity          | 90       |
| <b>OTHER FEATURES</b>              |                |                           |          |
| Size                               | Size 00        |                           |          |