

MY CARDIOMATE - THE POCKET ECG

Ref. 260080



CE
0197

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE

GARANTIE
2
2 YEARS
GUARANTEE



21/07/2026 00:06:43 – For information purposes only, can be subject to change

MY CARDIOMATE - THE POCKET ECG

Ref. 260080

Understand your heart in 30 seconds with a portable medical ECG.

Description

One click, 30 seconds, a complete ECG. Compact and easy to carry, My Cardiomate supports healthcare professionals in their daily practice by providing fast, simple, and reliable cardiac monitoring, wherever they are.

EASY TO USE

- Start the measurement with a single click
- Complete medical ECG in 30 seconds
- Built-in OLED screen that clearly displays your heart rate in real time and the ECG waveform.

ULTRA-PORTABLE DESIGN

- Fits in a pocket
- Weighs only 30 gr
- Very long battery life with 500 uses
- On-the-go use: heart monitoring anywhere.

COMPREHENSIVE DIGITAL ECOSYSTEM

- Free mobile app available on iOS and Android
- Advanced AI ECG analysis to identify abnormal ECG segments, AFib, and other arrhythmias in 30 seconds
- Unlimited data storage with cloud backup.
- Clinically approved in accordance with standard NF EN 60601-2-47 (safety and essential performance of ambulatory ECG systems)

MY CARDIOMATE - THE POCKET ECG

Ref. 260080

Features

CE Marking	CE0197	Medical device	Class IIa
Bar code	6934440102175	HS Code	9018110000
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	6	Packed unit / Height (cm)	17
Packed unit / Length (cm)	11	Packed unit / Weight	0.3
Bare Unit / Weight (kg)	0.04		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-25	Max temperature	70
Minimum humidity	10	Maximum humidity	95

OTHER FEATURES