


ELECTRONIC ARM BLOOD PRESSURE MONITOR NEOTENS® SPENGLER

Ref. NEO300



📅 27/07/2026 00:10:40 – For information purposes only, can be subject to change

ELECTRONIC ARM BLOOD PRESSURE MONITOR NEOTENS® SPENGLER



Ref. NEO300

NEOTENS® : ergonomic and precise arm blood pressure monitor, designed for daily use

Description

The **NEOTENS®** by Spengler is an **arm blood pressure monitor** that combines **ergonomics, precision, and patient comfort**. Its large LCD display (50 x 50 mm) ensures excellent readability, while the **upward measurement technology** delivers fast and comfortable readings at every use.

Thanks to its robust design and advanced functions, it is perfectly suited for healthcare professionals as well as patients requiring regular monitoring.

It is clinically approved for monitoring blood pressure during pregnancy.

Spengler's NEOTENS® is also recommended by STRIDE BP.

- **Display** : extra-large LCD screen (50 x 50 mm)
- **Measurement method** : Oscillometric, upward measurement for optimal comfort
- **Measurement site** : Arm
- **Range** :
Systolic pressure : 60 - 260 mm Hg
Diastolic pressure : 40 - 200 mm Hg
Pulse : 30 - 199 beats/min
- **Accuracy** : ± 3 mm Hg (pressure) ; ± 5 % (pulse)
- **Memory** : 50 readings per user x 2 users, with date and time stamp
- **Advanced features** :
Irregular heartbeat detection
Low battery indicator
Automatic shut-off after 1 minute of inactivity
- **Power supply** : 4 x AA batteries (≈ 300 measurements cycles) or optional USB-C power adapter

ELECTRONIC ARM BLOOD PRESSURE MONITOR NEOTENS® SPENGLER



Ref. NEO300

- **Cuff** : adult M/L size included (22-42 cm) ; child and XL adult cuffs sold separately
- **Robustness** : IP21 protecting rating
- **Weight & dimensions** : 230 g (without batteries), 108.6 x 103.5 x 73.5 mm
- **Medical classification** : Class IIa device (CE marking)
- **Warranty** : 3 years (excluding consumables)
- **Compatibility** : **SPENGLER HTA mobile app** (manual entry or photo) for simplified monitoring

The NEOTENS® offers **comfortable measurement, easy data tracking , and accuracy**, making it an essential tool for accurate and professional blood pressure monitoring.

Box Contents :

- **1 NEOTENS® arm blood pressure monitor**
- 1 M/L adult cuff (22-42 cm)
- 4 x AA alkaline batteries
- 1 carrying pouch
- 1 user manual

ELECTRONIC ARM BLOOD PRESSURE MONITOR NEOTENS® SPENGLER



Ref. NEO300

Features

CE Marking	CE	Medical device	Class IIa
Bar code	6942284608524	HS Code	90189010
Warranty (months)	36	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	11	Packed unit / Height (cm)	8
Packed unit / Length (cm)	28	Packed unit / Weight	0.855
Bare Unit / Weight (kg)	0.23		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-25	Max temperature	60
Minimum humidity	NA	Maximum humidity	93

OTHER FEATURES