

COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03_RC

ELITE BAGS



29/07/2026 23:46:35 – For information purposes only, can be subject to change

COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03_RC



Stylish kit with different compartments and elastic bands to carry the diabetic kit in an organized manner

Description

- 1 mesh bag (test strips, lancets, etc)- 1 pocket with velcro closure- 1 elastic for lancet device
- 2 elastics of different sizes (for pots, control solution, etc)
- 1 elastic with screen protector for the glucose reader- 1 insulated compartment for insulin pens with cold gel pocket
- 1 flat pocket outside
- Long shooters to facilitate opening
- Material: Polyester graybitone
- Weight: 0.12 kg.
- Measures: 17 x 10 x 6 cm / 0.95L
- Color: Grey (MEDIAB03), black (MEDIAB03_N), blue jeans (MEDIAB03_B) or pink (MEDIAB03_RC)
- Sold empty

COOL BAG FOR DIABETIQUE PM DIA, GREY

Ref. MEDIAB03_RC



Features

CE Marking	NONE	Medical device	No
Bar code	8436564542465	HS Code	4202.12.99.90
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	19	Bare Unit / Width (cm)	6
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	10
Packed unit / Length (cm)	10	Bare Unit / Length (cm)	17
Packed unit / Weight	0.248	Bare Unit / Weight (kg)	0.012
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	30	Maximum humidity	85
OTHER FEATURES			
Colors	Pink	Packed unit / Quantity	1