


MEDIC'S ELITE BAGS, BLUE

Ref. MEMEDI02_B

ELITE BAGS



01/08/2026 01:10:33 – For information purposes only, can be subject to change

MEDIC'S ELITE BAGS, BLUE

Ref. MEMEDI02_B



Description

Sanitary sports case with practical and comfortable design. Its internal organization totally adaptable thanks to the extractable velcro dividers, allows the rapid location of the equipment in case of medical or sports emergency.

- Quilted handle with velcro closure- Card holder
- Reinforced structure with ribbon
- Isothermal amplestick with 41 elastic strips- Document holder
- Large capacity external pocket with elastic bands and spacious compartments
- Adjustable honeycomb structure interior
- Flat pocket for the document holder
- 2 flat net pockets for instruments
- External pocket for personal oblets
- Measures: 50 x 25 x 31 cm
- Capacity: 33.6L
- Weight: 2.23 kg
- Material: 1200D polyester. Washable
- Recommended max. capacity: 12 kg
- Optional trolley
- Sold empty
- Available in blue (MEMEDI02_B) and black (MEMEDI02_N)

MEDIC'S ELITE BAGS, BLUERef. MEMEDI02_B**Features**

CE Marking	NONE	Medical device	No
Bar code	8436564540423	HS Code	4202290000
Warranty (months)	12	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	25
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	31
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	50
Packed unit / Weight	3.55	Bare Unit / Weight (kg)	2.23

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Colors	Blue	Packed unit / Quantity	1
--------	------	------------------------	---