

=====

Ref. E1.420.012

LUXAMED®



CE



 05/08/2026 22:27:07 – For information purposes only, can be subject to change

=====

Ref. E1.420.012

=====

LUXAMED®

Description

LUXAMED® Miller laryngoscope blades with integrated fiber optics are available in several sizes.

	Length (mm)	Width (mm)
Size 00	65	11
Size 0	77	11
Size 1	102	11
Size 2	155	14
Size 3	195	15
Size 4	210	18

- Blade for integrated fiber optics
- Made of high-quality stainless steel for a long service life
- Easy to clean and sterilize as there are no crevices to trap dust
- Maintenance-free
- 5,500 fibers for maximum service life and light intensity
- Optimal illumination with the largest possible cross-section (4 mm diameter)
- Illuminance of approximately 10,500 lux for 3.7V and 2.5V LEDs
- Compatible with all handles compliant with ISO 7376 (eco-friendly standard)
- Autoclavable up to 134°C (approximately 4,000 times)
- Data matrix with all important information for traceability

=====

Ref. E1.420.012

LUXAMED®

Features

CE Marking	CE	Medical device	Class I
Bar code	4260381960577	HS Code	90189084
Country	Germany		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	40
Minimum humidity	10	Maximum humidity	90
OTHER FEATURES			
Size	Size 0		