

## DIAGNOSTIC PENLIGHT - WITH STANDARD BULB - LUXAMED

Ref. D1.212.212

**LUXAMED®**



06/08/2026 22:53:24 – For information purposes only, can be subject to change

## DIAGNOSTIC PENLIGHT - WITH STANDARD BULB - LUXAMED

**LUXAMED®**Ref. D1.212.212

Penlight for eye and throat diagnostic., long life technology.

### Description

- Top quality ABS plastic housing
- Warm white light : comfortable for eyes
- Operating voltage 3.0 V
- Incl. 2 AAA batteries

Available in :

- White: D1.212.912
- Blue: D1.212.212
- Grey: D1.212.312
- Yellow: D1.212.512
- Kids : D1.312.912
- Black: D D1.212.112
- Red : D1.212.412
- Green : D1.212.612

**DIAGNOSTIC PENLIGHT - WITH STANDARD  
BULB - LUXAMED****LUXAMED®**Ref. D1.212.212**Features**

|                   |                |                |          |
|-------------------|----------------|----------------|----------|
| CE Marking        | CE             | Medical device | Class I  |
| Bar code          | 04260381961758 | HS Code        | 90189084 |
| Warranty (months) | 24             | Country        | Germany  |

**DIMENSIONS AND WEIGHT**

|                           |      |                           |    |
|---------------------------|------|---------------------------|----|
| Packed unit / Width (cm)  | nc   | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc   | Packed unit / Weight      | nc |
| Bare Unit / Weight (kg)   | 0.04 |                           |    |

**SPECIFIC STORAGE CONDITIONS**

|                     |     |                  |    |
|---------------------|-----|------------------|----|
| Minimum temperature | -10 | Max temperature  | 40 |
| Minimum humidity    | 10  | Maximum humidity | 90 |

**OTHER FEATURES**

|        |      |
|--------|------|
| Colors | Blue |
|--------|------|