

ISOTHERMAL BAG ELITE BAG MINICOOL'S, BLUE

Ref. MECOOL04

ELITE BAGS



13/08/2026 00:28:12 – For information purposes only, can be subject to change

ISOTHERMAL BAG ELITE BAG MINICOOL'S, BLUE

Ref. MECOOL04



Description

Insulated bag that combines the portability and practicality of ROW'S, with the capacity of COOL'S

The bag is a perfect size for laboratories transporting medium quantities of samples at home, as well as for samples from hospitalized patients, and even for vaccination campaigns.

- 1 isothermal main compartment with 4 terry supports to accommodate 44 tubes of different diameters, with 4 hermetic pockets to prevent liquid leakage.
- 1

insulated compartment for vials, with removable separators that prevent spilling•

The insulated compartments have pockets to place cold accumulators• 1

compartment with transparent front side for laboratory equipment•

1 external digital thermometer to monitor the temperature inside•

Outside, a transparent pocket for documents•

Stackable bag

- Material: 420D Polyester
- Outdoor measures: 45 x 24 x 27 cm / 29.2L•

Weight: 0.48 kg

- Recommended max. capacity: 12 kg

ISOTHERMAL BAG ELITE BAG MINICOOL'S, BLUE



Ref. MECOOL04

Features

CE Marking	NONE	Medical device	No
Bar code	3661040167459	HS Code	4202921900
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---