



SPO2 SENSOR ADULT FINGER CLIP FOR MASTER PALM 2

Ref. 222250



CE
0537



Accurate and reliable SpO₂ monitoring for adults, designed for use with the Spengler Master Palm 2.

Description

SpO₂ adult rigid sensor for MASTER PALM 2 oximeter

18/08/2026 13:34:04 – For information purposes only, can be subject to change

SPO2 SENSOR ADULT FINGER CLIP FOR MASTER PALM 2

Ref. 222250

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446030209	HS Code	90181990
Warranty (months)	12	Country	Uruguay
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.103
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	1	Max temperature	50
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		