

SPO2 SENSOR WRAP FOR MASTER PALM 2

Ref. 222253



CE
0537

DISPOSITIF MEDICAL
CLASSE
II_b
MEDICAL DEVICE

GARANTIE
1 AN
1 YEAR
GUARANTEE

SpO₂ wrap sensor for Master Palm 2 – for comfortable and reliable oxygen saturation monitoring.

Description

SpO₂ wrap sensor for oximeter MASTER PALM 2

18/08/2026 12:39:10 – For information purposes only, can be subject to change

SPO2 SENSOR WRAP FOR MASTER PALM 2

Ref. 222253

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446030230	HS Code	90181990
Warranty (months)	12	Country	Uruguay

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.085

SPECIFIC STORAGE CONDITIONS

Minimum temperature	1	Max temperature	50
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---