

SPO2 SENSOR NELLCOR OXIMAX COMPATIBLE NEONATE WRAP FOR MASTER PALM 3



Ref. 222263



CE
0537



18/08/2026 12:15:53 – For information purposes only, can be subject to change

SPO2 SENSOR NELLCOR OXIMAX COMPATIBLE NEONATE WRAP FOR MASTER PALM 3



Ref. 222263

Gentle and accurate oxygen monitoring for newborns, fully compatible with your Master Palm 3.

Description

NELCORR OXIMAX compatible SpO2 sensor wrap for oximeter MASTER PALM 3

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446030162	HS Code	90181990
Warranty (months)	12	Country	Uruguay

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.085

SPECIFIC STORAGE CONDITIONS

Minimum temperature	1	Max temperature	50
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---