


SPO2 SENSOR NELLCOR OXIMAX COMPATIBLE PEDIATRIC SILICON FOR MP3

Ref. 222265



CE
0537



18/08/2026 13:35:38 – For information purposes only, can be subject to change

SPO2 SENSOR NELLCOR OXIMAX COMPATIBLE PEDIATRIC SILICON FOR MP3



Ref. 222265

Description

NELCORR OXIMAX compatible pediatric silicon spo2 sensor for oximeter MASTER PALM 3

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3700446030186	HS Code	90181990
Warranty (months)	12	Country	Uruguay
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.084
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	1	Max temperature	50
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		