


MASTER PALM 3 OXYMETRE PORTABLE NELLCOR OXIMAX INSIDE

Ref. 223001



CE
0123



18/08/2026 13:33:11 – For information purposes only, can be subject to change

MASTER PALM 3 OXYMETRE PORTABLE NELLCOR OXIMAX INSIDE



Ref. 223001

The MASTER PALM 3® pulse oximeter was developed to meet the needs of healthcare professionals. Featuring Nellcoret NellcorOximax technology, it is condensed with technology.

Description

The Master Palm 3® pulse oximeter is a suitable material for healthcare professionals compatible with NELLCOR OXIMAX technology.

Its technology makes it possible to intervene quickly on all types of patients from adults to neonat, in order to accurately record their oxygen saturation and their frequency of hens.

Lightness and exceptional autonomy, you will appreciate the speed and precision of its measurements combined with its versatility.

• Adult and neonatal

technology • Exceptional autonomy • HD display with auto-rotate function that optimizes space utilization and readability by providing a vertical or horizontal layout.

• 2 display modes: Numbers and SpO2 Curve • Pulse, plethysmography curve, SpO2 and

SpO2 curve display • Fully configurable, it allows you to adjust the alarm, display, sound, system settings • Monitoring function (alarm and memories)

• Lightweight and compact 58.5 x 123 x 28 mm • Comes with rigid adult SpO2 sensor and 3 AA

batteries Accessories are available in addition according to your preferences: rechargeable battery and charging pad, pediatric or neonatal sensors, flexible sensors or clip, carrying case, etc ...

MASTER PALM 3 OXYMETRE PORTABLE NELLCOR OXIMAX INSIDE



Ref. 223001

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3700446027124	HS Code	9018199000
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.696
Bare Unit / Weight (kg)	0.105		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-20	Max temperature	55
Minimum humidity	15	Maximum humidity	85

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--