

SILICON PROTECT COVER FOR MASTER PALM 2/3 (2 PCS)

Ref. 223004



Flexible and reliable protection to extend the lifespan of your Master Palm 2 or 3 blood pressure monitor.

Description

Silicone protection for MASTER PALM 2 and MASTER PALM 3 oximeters.
Set of 2 protections, one yellow and one transparent.

18/08/2026 13:57:14 – For information purposes only, can be subject to change

SILICON PROTECT COVER FOR MASTER PALM 2/3 (2 PCS)

Ref. 223004



Features

CE Marking	NONE	Medical device	No
Bar code	3700446025434	HS Code	90330000
Country	China		

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.074

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---