

PAIR OF BLACK EARTIPS FOR PULSE II

Ref. 506900



18/08/2026 13:38:33 – For information purposes only, can be subject to change

PAIR OF BLACK EARTIPS FOR PULSE II

Ref. 506900



Description

Pair of standard-sized, black silicone ear tips for PULSE II and MAGISTER II

Features

CE Marking	NONE	Medical device	No
Bar code	3700446029982	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.004
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		