

## PAIR OF BLACK EARTIPS (SMALL) FOR PULSE II

Ref. 506901



18/08/2026 12:36:48 – For information purposes only, can be subject to change

## PAIR OF BLACK EARTIPS (SMALL) FOR PULSE II

Ref. 506901



Small ear tips designed for optimal compatibility with the Pulse II Magi stethoscope.

### Description

Pair of small, black silicone eartips for SPENGLER PULSE II and MAGISTER II

### Features

|                                    |               |                           |          |
|------------------------------------|---------------|---------------------------|----------|
| CE Marking                         | NONE          | Medical device            | No       |
| Bar code                           | 3700446029999 | HS Code                   | 90181990 |
| <b>DIMENSIONS AND WEIGHT</b>       |               |                           |          |
| Packed unit / Width (cm)           | nc            | Packed unit / Height (cm) | nc       |
| Packed unit / Length (cm)          | nc            | Packed unit / Weight      | 0.004    |
| <b>SPECIFIC STORAGE CONDITIONS</b> |               |                           |          |
| Minimum temperature                | NA            | Max temperature           | NA       |
| Minimum humidity                   | NA            | Maximum humidity          | NA       |
| <b>OTHER FEATURES</b>              |               |                           |          |
| Packed unit / Quantity             | 1             |                           |          |