


PAIR OF BLACK EARTIPS SILICONE FOR CARDIOPRESTIGE II

Ref. 506905



18/08/2026 13:36:32 – For information purposes only, can be subject to change

PAIR OF BLACK EARTIPS SILICONE FOR CARDIOPRESTIGE II

Ref. 506905



Soft and hygienic eartips for optimal comfort during every examination.

Description

Pair of SPENGLER CARDIO PRESTIGE II stethoscope ear tips.
 Black color, standard size, nano silicone.

Features

CE Marking	NONE	Medical device	No
Bar code	3700446030032	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.004
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		