

SPO2 SENSOR ADULT CLIP FOR COMEN NC3/NC5/C30

Ref. DBCADUL01



CE
0537

DISPOSITIF MEDICAL
CLASSE
II_b
MEDICAL DEVICE

GARANTIE
1 AN
1 YEAR
GUARANTEE

18/08/2026 12:17:50 – For information purposes only, can be subject to change

SPO2 SENSOR ADULT CLIP FOR COMEN NC3/NC5/C30

Ref. DBCADUL01

Designed for accurate and fast measurement of oxygen saturation

Description

Technical specifications: % oxygen: 90-100% accuracy +/-
-2% 70-90%

accuracy +/-3%
<70% not specified

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3661040219417	HS Code	90181990
Warranty (months)	12	Country	Uruguay

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	13	Packed unit / Height (cm)	3
Packed unit / Length (cm)	13	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---