


ADULT SPO2 SILICON SENSOR FOR COMEN NC3/NC5/C32

Ref. DBSADUL01

COMEN



CE
0537



ADULT SPO2 SILICON SENSOR FOR COMEN NC3/NC5/C32

COMEN

Ref. DBSADUL01

Designed to ensure accurate and fast measurement of oxygen saturation.

Description

Technical specifications: % oxygen: 90-100% Accuracy +/-2% 70-

90% Accuracy +/-3%
<70% not specified.

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3661040219431	HS Code	90181990
Warranty (months)	12	Country	Uruguay

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	13	Packed unit / Height (cm)	3
Packed unit / Length (cm)	13	Packed unit / Weight	nc

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---