

OXYMETER FINGER TYPE WITH CHILD SENSOR "DEPORTED" SOFT PC-60E

Ref. DPOXY02



CE
0123



📅 18/08/2026 12:37:50 – For information purposes only, can be subject to change

OXYMETER FINGER TYPE WITH CHILD SENSOR "DEPORTED" SOFT PC-60E

Ref. DPOXY02



For use in hospitals, doctors' offices and homes, for adults and children.

Description

The PC-60E suitable for use in the hospital, medical office and home.

- Automatic activation upon insertion of the finger into the oximeter
- Automatic shutdown after 6 seconds without activity
- OLED display resolution of SpO2, PR, bargraph and plethyspogramme
- Low battery indicator Visual
- Color Display Multi
- Audible and visual alarms
- Suitable for infants with Velcro sensor
- Dimensions: 66 x 36 x 33mm
- Weight: 60g (including batteries)
- Battery life: 30 hours continuous use
- Supplied with 2 AAA batteries, wrist strap 1 and 1 storage bag

OXYMETER FINGER TYPE WITH CHILD SENSOR "DEPORTED" SOFT PC-60E

Ref. DPOXY02



Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	6941900603110	HS Code	9018909911
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	15	Bare Unit / Width (cm)	3.1
Packed unit / Height (cm)	5	Bare Unit / Height (cm)	3.3
Packed unit / Length (cm)	23	Bare Unit / Length (cm)	6.1
Packed unit / Weight	0.1	Bare Unit / Weight (kg)	0.059
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		