

CHILD SENSOR "DEPORTED" SOFT FOR OXIMETER PC-60E

Ref. DPOXY04



CE
0123



18/08/2026 11:58:59 – For information purposes only, can be subject to change

CHILD SENSOR "DEPORTED" SOFT FOR OXIMETER PC-60E

Ref. DPOXY04



Soft pediatric SpO₂ sensor for PC-60E, designed for comfortable and reliable measurement in children

Description

With its remote pediatric sensor, flexible, reusable and immersible, it will adapt to children from 2 to 8 years old. The PC 60E pulse oximeter is reliable and accurate. It for health professionals in the city.

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040163833	HS Code	9018909911
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	14	Packed unit / Height (cm)	2
Packed unit / Length (cm)	14	Packed unit / Weight	0.001
Bare Unit / Weight (kg)	0.059		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--