

CHILD SENSOR "DEPORTED" VELCRO FOR OXIMETER PC-60E

Ref. DPOXY05



CE
0123



CHILD SENSOR "DEPORTED" VELCRO FOR OXIMETER PC-60E

Ref. DPOXY05



Soft neonatal SpO₂ sensor with VELCRO strap, designed for reliable measurement in infants.

Description

With its remote pediatric sensor, disposable VELCRO. It will adapt to the feet of infants from 0 to 2 years old. The PC 60E pulse oximeter is reliable and accurate. Ideal for urban healthcare professionals.

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040165998	HS Code	9018909911
Warranty (months)	6	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	15	Packed unit / Height (cm)	1
Packed unit / Length (cm)	14	Packed unit / Weight	0.059
Bare Unit / Weight (kg)	0.059		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--