

SET OF 5 SUNTECH MONOBLOC REUSABLE ARMBANDS ADULT PLUS (28-40CM) - BAYONET CONNECTOR

SUNTECH

Ref. DSU231016A



CE



Description

Set of 5 Suntech adult plus nylon armbands (28-40cm)

18/08/2026 13:56:08 – For information purposes only, can be subject to change

SET OF 5 SUNTECH MONOBLOC REUSABLE ARMBANDS ADULT PLUS (28-40CM) - BAYONET CONNECTOR

SUNTECHRef. DSU231016A

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446013233	HS Code	9018905040
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	23	Packed unit / Height (cm)	9
Packed unit / Length (cm)	25	Packed unit / Weight	0.584

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---