

## FLEXIBLE CHILD DIGITAL SPO2 SENSOR (90CM CABLE)

**SUNTECH**Ref. DSU250904

Reliable and gentle oxygen monitoring for children, compatible with SunTech devices.

### Description

Flexible child digital SpO2 sensor (90 cm cable)

### Features

|            |               |                |           |
|------------|---------------|----------------|-----------|
| CE Marking | CE0413        | Medical device | Class IIb |
| Bar code   | 8698870322296 | HS Code        | 90181990  |
| Country    | Turkey        |                |           |

#### DIMENSIONS AND WEIGHT

|                           |    |                           |      |
|---------------------------|----|---------------------------|------|
| Packed unit / Width (cm)  | 2  | Packed unit / Height (cm) | 16   |
| Packed unit / Length (cm) | 14 | Packed unit / Weight      | 0.01 |

#### SPECIFIC STORAGE CONDITIONS

|                     |    |                  |    |
|---------------------|----|------------------|----|
| Minimum temperature | NA | Max temperature  | NA |
| Minimum humidity    | NA | Maximum humidity | NA |

#### OTHER FEATURES

|                        |   |
|------------------------|---|
| Packed unit / Quantity | 1 |
|------------------------|---|