

DISPOSABLE PROBE COVERS FOR THERMOSCAN PRO 6000, BOX OF 800

Ref. DWPRO6K02

WelchAllyn®

Advancing Frontline Care™



CE



18/08/2026 12:15:26 – For information purposes only, can be subject to change

DISPOSABLE PROBE COVERS FOR THERMOSCAN PRO 6000, BOX OF 800

Ref. DWPRO6K02

WelchAllyn®
Advancing Frontline Care™

Disposable probe covers for ThermoScan Pro 6000, ensuring hygiene and reliability with every reading

Description

This covers probe helps protect your ThermoScan Pro 6000 ear thermometer.

Single use for better efficiency and impeccable hygiene, this protective cap is perfectly suited to the Braun thermometer.

Its shape is studied to suit children and babies.

Colour: transparent

Features

CE Marking	CE	Medical device	Class I
Bar code	00732094321043	HS Code	9025900090
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	8
Packed unit / Length (cm)	30	Packed unit / Weight	0.5
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		