

THERMOSCAN PRO 6000 RECHARGEABLE BATTERY PACK 1

Ref. DWPRO6K09

WelchAllyn®

Advancing Frontline Care™



Description

- 2.4V 750mAh rechargeable battery for Pro 6000 thermometer
- Fully charged, the thermometer can make 700 measurements
- Must be recharged thanks to the charging dock

📅 18/08/2026 13:38:56 – For information purposes only, can be subject to change

THERMOSCAN PRO 6000 RECHARGEABLE BATTERY PACK 1

Ref. DWPRO6K09

WelchAllyn®

Advancing Frontline Care™

Features

CE Marking	CE	Medical device	No
Bar code	00732094325089	HS Code	850610
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	2	Packed unit / Height (cm)	1
Packed unit / Length (cm)	4	Packed unit / Weight	0.042
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		