

OXYSTART® - DIGITAL FINGERTIP PULSE OXIMETER BLUBERRY

Ref. 221052



CE
0123



18/08/2026 13:31:25 – For information purposes only, can be subject to change

OXYSTART® - DIGITAL FINGERTIP PULSE OXIMETER BLUBERRY



Ref. 221052

For SpO2 spot check.

Description

HIGH MEASUREMENT

ACCURACY• Plethysmography curve display•
High sensitivity sensor

EASY TO USE AND

CONVENIENT• Compact and lightweight•
Switch on at the touch of the clock•
Automatic shutdown after 5 seconds of non-use•
Automatic display
rotation• Simultaneous display in less than 10 sec. :
-pulsed oxygen saturation (% SpO2),
- heart rate,
- the intensity of the pulsatile signal (graph bar),
- the plethysmography curve.
• TFT display•
Low battery
indicator• Customizable display: 6 display modes.
• 5 levels of brightness

ADDITIONAL

INFORMATION• Comes with a wrist strap, a cover and 2 AAA LR03 batteries

OXYSTART® - DIGITAL FINGERTIP PULSE OXIMETER BLUBERRY



Ref. 221052

Features

CE Marking	CE0123	Medical device	Class IIb
Bar code	3700446041809	HS Code	90181910
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.135
Bare Unit / Weight (kg)	0.06		

SPECIFIC STORAGE CONDITIONS

Minimum temperature	-40	Max temperature	60
Minimum humidity	5	Maximum humidity	95

OTHER FEATURES

Packed unit / Quantity	1		
------------------------	---	--	--